BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 1 of 16	
Colour Reference:	BLACK	

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Date of Revision: 2019/12/01

Baxter

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 2 of 16	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 3 of 16	
Colour Reference:	BLACK	

◆ eng	English	4
◆ chi	繁體中文	6
◆ kor	한국어	8
	General characteristics.....	11
	Installation guide	14

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 4 of 16	
Colour Reference:	BLACK	

ENGLISH

eng

The **PRISMAFLEX HF20 Set** is manufactured by **GAMBRO Industries**,
7 Avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

⚠ Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

In this document:

⚠ "Warning" is used to alert the user/operator not to take a certain action which, if taken, can cause a potential hazard and result in a serious adverse reaction, injury or death.

⚠ "Caution" is used to alert the user/operator to take a certain action to protect against a potential hazard, which, if ignored, could have an adverse effect on the patient or on the device.

"Note" is used as a reminder to the user/operator on normal treatment activity and on what is a suitable action in a particular situation.

SCUF: Slow Continuous UltraFiltration.

CVVH: Continuous Venous Hemofiltration.

CVVHD: Continuous Venous HemoDialysis.

CVVHDF: Continuous Venous HemoDiaFiltration.

Prdilution: addition of replacement fluid to the blood stream upstream to the filter.

Postdilution: addition of replacement fluid to the blood stream downstream to the filter.

PRODUCT DESCRIPTION

- The **PRISMAFLEX HF20 Set** is a disposable, extracorporeal circuit for use with the **PRISMAFLEX** system equipped with software 4.00 or later.
- The **PRISMAFLEX HF20 Set** consists of a PAES hollow fiber hemofilter/dialyzer* and tubing lines; refer to **PRISMAFLEX** control unit operator manual drawings for details.
- All line connectors are compatible with the ISO 594/1 & 2 international standards concerning conical fittings.
- The fluid pathways of the **PRISMAFLEX** Set are guaranteed sterile and non pyrogenic.
- The **PRISMAFLEX HF20 Set** is sterilized by ethylene oxide (EtO). Deaeration is such that EtO residuals comply with the ones described in ISO 10993-7.
- Expiration date: Please refer to product label.

* In this document the hemofilter/dialyzer will be referred to as "filter".

INTENDED USE / INDICATIONS

The **PRISMAFLEX** Set is indicated for use only with the **PRISMAFLEX** control unit equipped with software version 4.00 or later in providing continuous fluid management and renal replacement therapies. The system is intended for patients who have acute renal failure, fluid overload, or both.

This set is intended for use in the following veno-venous therapies: SCUF; CVVH; CVVHD; CVVHDF.

All treatments administered via the **PRISMAFLEX** Set must be prescribed by a physician.

The size, weight, state of uremia, cardiac status, and general physical condition of the patient must be carefully evaluated by the prescribing physician before each treatment.

CONTRAINDICATIONS

There are no known absolute contraindications for the **PRISMAFLEX HF20** device. For the following conditions a careful assessment of the individual patient's risk/benefit ratio has to be made by the treating physician (relative contraindications): Inability to establish vascular access; Inability to provide adequate or safe anticoagulation (e.g., uncontrolled bleeding diathesis, central nervous system bleed); Severe hemodynamic instability.

CAUTIONS AND WARNINGS

Note: additional warnings and cautions pertaining to the **PRISMAFLEX** system are included in the **PRISMAFLEX** control unit operator's manual.

⚠ Cautions

It is recommended that particular attention must be paid with respect to extracorporeal blood volume; as the 'extracorporeal volume/patient blood volume' ratio is very high, the physician may decide to prime the extracorporeal circuit with adequate volume substitution before patient connection. Refer to the "Before connecting patient..." section of this document and to General characteristics.

The **PRISMAFLEX HF20 Set** should be restricted to patients with a body weight greater than 8kg (18lb).

- Carefully read these instructions for use and the **PRISMAFLEX** control unit operator's manual before using this product.
- Store the **PRISMAFLEX** Set in a dry place, between 0° C (32°F) and 30°C (86°F).
- Some solvents and other chemicals, if used in contact with the filter, could damage the set. No chemical of this type should be used without permission of the manufacturer. The following are especially forbidden:
 - halogenated aromatic and aliphatic solvents,
 - ketonic solvents.
- To prevent contamination, this **PRISMAFLEX** Set must be used as soon as its packaging and sterilization caps are removed.
- Do not use this set if the packaging is damaged, if the sterilization caps are missing or loose, or if any of the lines in the set are kinked.
- Do not try to remove the filter from the cartridge plate.
- Destroy this set after single use, using aseptic technique for potentially contaminated equipment. Do not resterilize. The **PRISMAFLEX HF20 Set** is intended for single use only. Re-using the **PRISMAFLEX HF20 Set** may cause serious damage to the product resulting in patient injury or death.
- Use aseptic techniques when handling all blood and fluid lines in the set.
- Connect the **PRISMAFLEX** Set to a patient via venous blood access and return devices. A double-lumen venous catheter is the recommended blood access device; however, two single-lumen venous catheters can also be used. There are 3 possible accesses for **PRISMAFLEX** System therapies: subclavian, jugular or femoral vein.
- If the patient is not immediately connected to the **PRISMAFLEX** Set after priming is complete, flush the set with at least 500 mL priming solution [saline or alkaline solution (pH ≥ 7.3) with heparin added] prior to connecting the patient. This requires use of a new bag of priming solution.
- Use a 21-gauge or smaller needle to obtain blood/fluid samples or remove trapped air from the **PRISMAFLEX** Set. Use of larger needles can cause holes in the sample sites, resulting in blood loss or air embolism.
- In the case of patients who pose a high risk of haemorrhaging, it is recommended not to add heparin to the priming solution.
- When not using the pre blood pump infusion circuit, it is recommended to clamp this circuit close to its connection to the access line; this will prevent the sedimentation of blood into the pre blood infusion line.
- When treating low weight patients (i.e. patient with a weight between 8 and 20 kgs) with the **PRISMAFLEX HF20 Set**, the volume of the extracorporeal circuit and the high 'exchange volume/patient weight' ratio may cause hypothermia. We therefore recommend warming the blood.
- The **PRISMAFLEX HF20 Set** is compatible with the sleeve blood warmers, which should be installed on the return line. Refer to the specific instructions for use.
- If citrated blood is used for the priming of the extracorporeal circuit before patient connection, it is recommended that the pH of this blood be verified and buffered at a value between 7.3 and 7.5. The ionized calcium should be adjusted at a value > or = 1.0.

⚠ Warnings

- The use of operating procedures other than those published by the manufacturer or the use of accessory devices not recommended by the manufacturer can result in patient injury or death.
- Should acute allergic reactions (first-use syndrome) occur in patients receiving treatment via the **PRISMAFLEX HF20 Set**, immediately stop the treatment and administer appropriate intervention.
- Do not allow air to enter the blood compartment of the filter after priming is started. If a large amount of air enters, the set must be replaced.
- To assure adequate filter performance, it is recommended that the set be changed every 24 hours of use. However, the set must be changed after 3 days (72 hours). Continued use beyond this limit could result in rupture of the pump segments, with risk of patient injury or death.
- Use the **PRISMAFLEX HF20 Sets** only with the **PRISMAFLEX** control unit equipped with software version 4.00 or later. The use of previous software versions can result in patient injury or death.

SPECIFICATIONS

See Tables at the end of this IFU.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 13-JAN-2020	Proofread No.: 2
Designer: CDS	Page: 5 of 16	
Colour Reference:	BLACK	

SET MATERIALS

PAES hollow fiber	: PolyaryletherSulfone
Housing and headers	: Polycarbonate
Potting compound	: Polyurethane
Tubing material	: Plasticized polyvinyl chloride (PVC)
Cartridge	: PETG

Note: the following information is available from the manufacturer upon request:

- information about test methods used to obtain performance characteristics,
- the types and amounts of residue from the sterilization process.

The **PRISMAFLEX** HF20 Set is not made with natural rubber latex.

INSTRUCTIONS FOR USE

Note: use the set by following the detailed on-line instructions provided by the **PRISMAFLEX** control unit. Additional information is available in the **PRISMAFLEX** control unit operator's manual.

Load Set

Install the set onto the **PRISMAFLEX** control unit using the photographs on the inside cover as a guide.

Prepare and Connect Solutions

Hang bag of priming solution [saline or alkaline solution (pH ≥ 7.3) with added 5000IU heparin/liter] on priming hook (left corner hook top of front panel). Connect access (red)/effluent (yellow) Y-line to priming solution bag.

Before connecting patient

If priming of circuit with adequate volume substitution is requested, refer to the **PRISMAFLEX** control unit's user interface for specific instructions.

Anticoagulation considerations

According to the literature* continuous heparinization at a rate of 10 to 20 IU/kg/h ensures proper operation of the extracorporeal circuit when performing treatment with patients having normal coagulation status.

Depending on the patient's condition, however, heparinization can be lowered to less than 5 IU/kg/h.

Heparinization can be controlled, for instance, by partial thromboplastin time (PTT) measurements: in this case, PTT could be maintained at 20 to 30 seconds over baseline.

SPECIAL PROCEDURES IN CASE OF COMPLICATION

External Blood Leaks

External blood leakage may not be immediately identified by monitoring equipment and could result in significant blood loss. Check the filter and all connections of the disposable tubings during treatment to minimize the risk of leakage.

* "Continuous renal replacement therapy in critically ill neonates", G.ZOBEL & Al., *Kidney International*, Vol.53, Suppl. 66 (1998).

If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the **PRISMAFLEX** Set.

If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

Hypersensitivity Reactions

Note: see Warning no. 2.

Should acute allergic reactions (first use syndrome) occur within the first few minutes of the treatment, it is important to react immediately by discontinuing the session and administering the appropriate treatment.

Adverse reactions may occur due to the complex interaction between blood and the artificial surfaces of the entire extracorporeal circuit. These reactions may also be precipitated and/or exacerbated by other external factors involved with the individual patient's specific disease process and the treatment of renal insufficiency. Certain types of adverse reaction may occur due to operational factors associated with the treatment. Therefore, proper management of the fluid removal, electrolyte balance, anticoagulation and blood flow rate as well as monitoring of the overall treatment parameters are essential to avoid side-effects which may be associated with hemodialysis/hemofiltration therapies.

Hypersensitivity reactions have been observed during dialysis. Symptoms of a hypersensitivity reaction may be gastrointestinal, mucocutaneous, respiratory, cardiovascular or systemic in nature and range from very mild to severe. Such symptoms have been described as anaphylactic-like reactions within the first few minutes. Manifestations include nausea, malaise, weakness, a sensation of burning or heat throughout the body, profuse perspiration, respiratory distress and in some instances hypotension and cardiopulmonary arrest. Should a combination of such symptoms appear, particularly at the start of the treatment session, it is important to react immediately by discounting the session and administering appropriate treatment. Blood in the extracorporeal circuit must not be returned to the patient.

Extra care must be taken when treating patients who have exhibited possible hypersensitivity symptoms during previous treatments, or patients who have a history of being highly sensitive and allergic to a variety of substances. A physician must be consulted to evaluate the risk and prescribe the appropriate precautions if a possible sensitivity is suspected.

The following factors are considered essential to minimize the risk of hypersensitivity reaction and other side effects:

1. Strict adherence to the set-up, priming and rinsing procedures detailed in the manufacturer's Instructions for Use.
2. Setting up and monitoring the treatment operating parameters according to the manufacturer's recommendations specified for each type of **PRISMAFLEX** Set and to the patient's needs and tolerance.
3. Strict adherence to all WARNINGS and CAUTIONS given by the manufacturer in the instructions for use.

WARRANTY AND LIMITATION OF LIABILITY

- a) The manufacturer warrants that the **PRISMAFLEX** Set has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.
If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the **PRISMAFLEX** Set becoming apparent before the expiration date.
- b) The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the **PRISMAFLEX** Set and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the **PRISMAFLEX** Set, whether as a result of any defect therein or otherwise.
- c) The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer's release of the **PRISMAFLEX** Set, failure or omission to inspect the **PRISMAFLEX** Set before use in order to ensure that the **PRISMAFLEX** Set is in proper condition, or any warranty given by independent distributors or dealers.
- d) The manufacturer is GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

INDIA

- **Name and address of importer:** Baxter India Pvt. Ltd., Gala No. 1 to 6, Building No. C15, Ground Floor, Shree Arihant complex, Kalher, Bhiwandi, Dist. Thane,
Tal: Bhiwandi (Thane-Zone 5) 421302
- Name, address, telephone no., email address (in case of consumer complaint):
Baxter India Pvt. Ltd.,
5th Floor, Block A, Building 9, DLF Phase III, Cyber City, Gurgaon- 122002.
- Consumer Care No.: 0124-4603200
- Consumer Care email ID: Customerservice_SHS_India@baxter.com
- Import License No.: MD-1591-2310

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 6 of 16	
Colour Reference:	BLACK	

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- SUCF (Slow Continuous Ultrafiltration)
- CVVH (Continuous Venovenous Hemofiltration)
- CVVHD (Continuous Venovenous Hemodialysis)
- CVVHDF (Continuous Venovenous Hemodiafiltration)

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
BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 10 of 16	
Colour Reference:	BLACK	

"In vitro" ultrafiltration with blood (values \pm 20 %).
 (Bovine blood at 37°C, Hct 32%, Protein concentration 60g/L).
 Ultrafiltration is controlled by the PRISMAFLEX system and is independent of the ultrafiltration coefficient (KUF).

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 11 of 16	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY

Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 12 of 16	
Colour Reference: BLACK		

FILTER DATA	過濾器資料	필터 데이터	PRISMAFLEX HF 20 SET
NOMINAL PHYSICAL CHARACTERISTICS Effective surface area Fiber internal diameter (wet) Fiber wall thickness	標稱物理規格 有效表面面積 纖維內徑 (潮濕) 纖維壁厚	물리적 공칭 특성 유효 표면적 섬유 내부 직경 (젖은 경우) 섬유벽 두께	0.2 m ² 215 μm 50 μm
IN VITRO PERFORMANCES ≠ Blood priming volume	體外效能 ≠ 血液預充量	시험관내 성능 ≠ 혈액 프라이밍량	17 mL
Blood pressure drop (post dilution) (bovine blood, Htc****32%, Cp***** 60 g/L, 37°C) QB** = 20 mL/min QB** = 50 mL/min QB** = 100 mL/min	血壓下降 (後稀釋) (牛血, Htc****32%, Cp***** 60 g/L, 37°C) QB** = 20 毫升/分鐘 QB** = 50 毫升/分鐘 QB** = 100 毫升/分鐘	혈액 압력 강하(희석 후) (소 혈액, Htc****32%, Cp***** 60g/L, 37°C) QB** = 20 mL/min QB** = 50 mL/min QB** = 100 mL/min	QUF*** = 0 L/h 25 mmHg 37 mmHg 59 mmHg
Sieving coefficient (bovine plasma, Cp 60 g/L, 37° C) QB** = 50 mL/min at 37°, QUF*** = 10 mL/min - Urea - Vitamin B12 - Inulin - Albumin	篩濾係數 (牛血漿, Cp 60 g/L, 37° C) QB** = 37° 時為 50 毫升/分鐘, QUF*** = 10 毫升/分鐘 - 尿素 - 維他命 B12 - 菊糖 - 白蛋白	Sieving 계수 (소 혈장, Cp 60g/L, 37°C) QB** = 50 mL/min (37°C), QUF*** = 10 mL/min - 요소 - 비타민 B12 - 이눌린 - 알부민	1.00 1.00 0.92 < 0.01






≠ Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions. / 一般由實驗室測試滅菌後樣品所得之平均值。依據病人及臨床狀況，結果可能不同。 /

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


BAXTER CONFIDENTIAL - INTERNAL USE ONLY	
Part Number: 1000014798	Date: 08-JAN-2020 Proofread No.: 1
Designer: CDS	Page: 13 of 16
Colour Reference:	BLACK

THIS PAGE HAS BEEN LEFT DELIBERATELY BLANK

INSTALL SET DIRECTLY ONTO **Prismaflex** **MONITOR**

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014798	Date: 08-JAN-2020	Proofread No.: 1
Designer: CDS	Page: 15 of 16	
Colour Reference:	BLACK	

Patented ready to use PrismaFlex set packaging.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY	
Part Number: 1000014798	Date: 08-JAN-2020 Proofread No.: 1
Designer: CDS	Page: 16 of 16
Colour Reference:	BLACK