BAXTER CONFIDENTIAL - INTERNAL USE ONLY	
Part Number: 1000014865	Date: 06-JAN-20 Proofread No.: P03
Designer: E.M.	Page: 1 of 24
Colour Reference:	BLACK

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
BAXTER CONFIDENTIAL - INTERNAL USE ONLY	
Part Number: 1000014865	Date: 06-JAN-20 Proofread No.: P03
Designer: E.M.	Page: 3 of 24
Colour Reference:	BLACK

English.....

India.....

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Installation guide

6

9

10

14

20

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 4 of 24	
Colour Reference:	BLACK	

DEFINITION OF SYMBOLS

		ENGLISH		繁體中文	
SYMBOL GRAPHIC & REF NUM / 符號圖形和參考編號	TITLE AND NUMBER OF STANDARD / 標題和標準編號	SYMBOL TITLE	SYMBOL DESCRIPTION (EXPLANATORY TEXT)	符號名稱	符號描述（解釋性文字）
5.1.1 	ISO 15223-1 ¹	Manufacturer	Indicates the medical device manufacturer, i.e. the natural or legal person with responsibility for the design, manufacture, packaging and labelling of a device before it is placed on the market under his own name.	製造廠	表示醫療器材製造廠，即在以其名義將器材上市前負責器材設計、生產、包裝和標籤的自然人或法人。
5.1.3 	ISO 15223-1 ¹	Date of manufacture	Indicates the date when the medical device was manufactured. Format should be YYYY-MM-DD.	製造日期	表示生產該醫療器材的日期。格式應為年-月-日。
5.1.4 	ISO 15223-1 ¹	Use-by date	Indicates the date after which the medical device is not to be used. Format should be YYYY-MM-DD	保存期限	表示在該日期後不應使用該醫療器材。格式應為年-月-日
5.1.5 	ISO 15223-1 ¹	Batch code	Indicates the manufacturer's batch code. Synonyms are "lot number" and "batch number"	批次代碼	表示製造廠的批次代碼。同義字為“批號”和“批次號”
5.1.6 	ISO 15223-1 ¹	Catalogue number	Indicates the manufacturer's catalog number. Synonyms are "reference number" and "reorder number"	目錄編號	表示製造廠的目錄編號。同義字為“參考編號”和“再訂購編號”
5.2.3 	ISO 15223-1 ¹	Sterile fluid path	Indicates the presence of a sterile fluid path within the medical device in cases where other parts of the medical device, including the exterior, might not be supplied sterile. The method of sterilization shall be indicated in the empty box, as appropriate.	無菌液體通路	表示當該醫療器材其它部分（包括外部）不能以無菌形式供應情況下，該醫療器材的內部為無菌液體通路。應將滅菌方法標示在外盒。
5.3.1 	ISO 15223-1 ¹	Fragile, handle with care	Indicates a medical device that can be broken or damaged if not handled carefully.	易碎、小心輕放	表示如果不注意操作，可能會打破或損壞醫療器材。
5.3.4 	ISO 15223-1 ¹	Keep dry OR Keep away from rain	Indicates a medical device that needs to be protected from moisture. OR Indicates a medical device that needs to be kept away from rain.	保持乾燥或遠離雨水	表示醫療器材需要防潮。或表示醫療器材需要遠離雨水。
5.3.7 	ISO 15223-1 ¹	Temperature Limit	Indicates the temperature limits to which the medical device can be safely exposed.	溫度限值	表示該醫療器材可安全地暴露於該溫度限值內。
5.4.2 	ISO 15223-1 ¹	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure. Synonyms for "Do no re-use" are "single use" and "use only once".	請勿重複使用	表示醫療器材僅供單次使用，或用於單一病人的單次治療。“請勿重複使用”的同義字為“單次使用”和“僅使用一次”。
5.4.4 	ISO 15223-1 ¹	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.	注意	表示使用者需要參閱使用說明、了解重要警訊，如因為各種原因不能標示在醫療器材機體上的各種警告和預防措施。
5.4.3 	ISO 15223-1 ¹	Consult instructions for use	Indicates the need for the user to consult the instructions for use. Synonym for "Consult instructions for use" is "Consult operating instructions".	參閱使用說明	表示使用者需要參閱使用說明。“參閱使用說明”的同義字為“參閱操作手冊”。
	NA, European countries only / 不適用，僅限歐洲國家	Green Dot Symbol	To indicate the manufacturer of the product contributes to the cost of recovery and recycling	綠點符號	表示該產品製造廠落實回收和重複利用
5.2.8 	ISO 15223-1 ¹	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or open	如果包裝破損，請勿使用	表示如果包裝已破損或打開，不應使用醫療器材
0623 	ISO 7000 ²	This way up	Indicates the correct, upright position of the package	此面朝上	表示包裝的正確、朝上位置

¹ ISO 15223-1:2016 Medical Devices – Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General Requirements / 醫療器材 - 要提供與醫療器材標籤配合使用的符號、標記和資訊第 1 部分：基本要求

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 6 of 24	
Colour Reference:	BLACK	

ENGLISH

eng

The PRISMAFLEX TPE1000/TPE2000 Set is manufactured by GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

DEFINITION OF EXPRESSIONS USED IN THIS MANUAL

In this document :

- ⚠ “**Warning**” is used to alert the user/operator **not to take** a certain action which, if taken, can cause a potential hazard and result in a serious adverse reaction, injury or death.
- ⚠ “**Caution**” is used to alert the user/operator **to take** a certain action to protect against a potential hazard, which, if ignored, could have an adverse effect on the patient or on the device.
- “**Note**” is used as a reminder to the user/operator on normal treatment activity and on what is a suitable action in a particular situation.

TPE Therapeutic Plasma Exchange.

PRODUCT DESCRIPTION

- The PRISMAFLEX TPE1000/TPE2000 Set is a disposable, extracorporeal circuit for use with the PRISMAFLEX System having selected ‘TPE’ therapy.
- The PRISMAFLEX TPE1000/TPE2000 Set consists of a polypropylene hollow fiber plasmafilter* and tubing lines.
- This filter is permanently connected to a blood access line (red-striped), a blood return line (blue-striped) and an effluent outlet line (yellow-striped).
- The other lines of the set include :
 - a replacement solution line (purple-striped),
 - a pre blood pump line (white striped),
 - an anticoagulant line (syringe).
- The pre blood pump line allows the addition of infusion solution close to the end of the patient access line and before the blood pump. This can be used as an additional pre-dilution infusion to the replacement circuit.
- The PRISMAFLEX Set is provided with a specific small volume deaeration chamber in which blood does not appear to mix with the replacement liquid the majority of the time; this is a normal operation of the device.
- A 5-liter bag is provided to be connected to the Y line connecting ends of the effluent and blood return lines to initially collect priming solution, during priming. Then, during treatment, this bag is used to collect filtrated plasma. Other sterile 5 and 9 liter bags and sterile, non pyrogenic spikes can be ordered separately.
- All line connectors are compatible with the ISO 594/1 & 2 international standards concerning conical fittings.
- The fluid pathways of the PRISMAFLEX Set are guaranteed sterile and non pyrogenic.
- The PRISMAFLEX TPE1000/TPE2000 Set is sterilized by ethylene oxide (EtO). Deaeration is such that EtO residuals comply with the ones described in ISO 10993-7.
- The PRISMAFLEX TPE1000/TPE2000 Set complies with ISO 13960.
- Expiration date: please refer to product label.

* In this document the plasmafilter will be referred to as “filter”.

INTENDED USE / INDICATIONS

The PRISMAFLEX TPE1000/TPE2000 Set is indicated for use only with the PRISMAFLEX control unit, after having selected TPE therapy (therapeutic plasma exchange).

The PRISMAFLEX TPE1000/TPE2000 Set is intended for use in therapeutic plasma exchange, thus in diseases where removal of plasma components is indicated.

The use of the PRISMAFLEX TPE1000 Set should be restricted to patients with a body weight greater than 9kgs (20lbs).

The use of the PRISMAFLEX TPE2000 Set should be restricted to adults.

The size, weight, state of uremia, cardiac status and general physical condition of the patient must be evaluated by the prescribing physician before each treatment.

The device should be used only on the direction of a physician who has evaluated all of the pertinent features of this device in relation to the individual patient.

Patients with bleeding tendencies must be closely supervised during the treatment. Treatment should be discontinued if acute hemorrhage that can not be corrected (hematemesis, hemoptysis and melena) occurs during TPE).

All treatments administered via the PRISMAFLEX Set must be prescribed by a physician.

CONTRAINDICATIONS

There are no known contraindications to therapeutic plasma exchange.

CAUTIONS AND WARNINGS

Note: additional warnings and cautions pertaining to the PRISMAFLEX system are included in the PRISMAFLEX operator’s manual.

⚠ Cautions

It is recommended that particular attention must be paid with respect to extra-corporal blood volume (see General Characteristics).

- Carefully read these instructions for use and the PRISMAFLEX operator’s manual before using this product.
- Store the PRISMAFLEX Set in a dry place, between 0° C (32° F) and 30° C (86° F).
- Some solvents and other chemicals, if used in contact with the filter, could damage the set. No chemical of this type should be used without permission of the manufacturer. The following are especially forbidden :
 - a) halogenated aromatic and aliphatic solvents,
 - b) ketonic solvents.
- To prevent contamination, this PRISMAFLEX Set must be used as soon as its packaging and sterilization caps are removed.
- Do not use this set if the packaging is damaged, if the sterilization caps are missing or loose, or if any of the lines in the set are kinked.
- Do not try to remove the filter from the cartridge plate.
- Destroy this set after single use, using aseptic technique for potentially contaminated equipment. Do not resterilize. The PRISMAFLEX TPE 1000/ TPE 2000 Set is intended for single use only. Re-using the PRISMAFLEX TPE 1000/ TPE 2000 Set may cause serious damage to the product resulting in patient injury or death.
- Use aseptic techniques when handling all blood and fluid lines in the set.
- Use only prescribed replacement solutions with the PRISMAFLEX system. **The prescribed solution must take into account the protein loss by plasmafiltration.**
The use of a nonsterile replacement solution could induce risks of bacterial or pyrogenic contamination for the patient.
- Connect the PRISMAFLEX Set to a patient via venous blood access and return devices. A double-lumen venous catheter is the recommended blood access device; however, two single-lumen venous catheters can also be used. There are 3 possible accesses for PRISMAFLEX system therapies : subclavian, jugular or femoral vein. If the patient is provided with an arterio-venous fistula access, the therapeutic plasma exchange treatment can be performed with standard dialysis needles.
- During priming and operation, observe closely for leakage at joints and connections within the set, notably the bags. Leakage can cause blood loss or air embolism. If leakage cannot be stopped by tightening the connections, replace the set.
- Before connecting the blood return line to the patient, check for absence of air between the segment of line inserted in the air detector and the patient-end of the return line.
If air is present in this part of the return line, connect the access line to the patient and start the blood pump while leaving the return line connected to the effluent bag. Purge the air present in the end-part of the return line, then stop the blood pump. Disconnect the return line from the effluent bag, and connect it to the patient.
If the amount of air in the blood circuit is too large, reprime the circuit completely before patient connection.
- After priming is complete, do not remove the pressure pods from the pressure sensor housings. If pods are removed, the set must be changed or the Diaphragm Reposition procedure performed (refer to PRISMAFLEX operator’s manual).
- If the patient is not immediately connected to the PRISMAFLEX Set after priming is complete, flush the set with at least 1 000 mL priming solution [saline or alkaline solution (pH ≥ 7.3) with heparin added] prior to connecting the patient. This requires use of a new bag of priming solution.
- Use a 21-gauge or smaller needle to obtain blood/fluid samples or remove trapped air from the PRISMAFLEX Set. Use of larger needles can cause holes in the sample sites, resulting in blood loss or air embolism.
- The PRISMAFLEX system may not be able to detect disconnections of the set from the patient’s catheter. Carefully observe the set and all operations while using the PRISMAFLEX system for a patient treatment.
- 17a.** Due to the nature of use of the PRISMAFLEX Set (low blood flow rates and other special factors), the possibility for coagulation within the blood flowpath is substantially enhanced. Pay careful attention to the possible medical hazards associated with coagulation of the blood flowpath and comply with the minimum blood flow rates specifications of each filter (see the “Filter Operating Specifications” section).
- 17b.** Reference literature* and laboratory data have shown an increase in the

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 7 of 24	
Colour Reference:	BLACK	

plasmafilter pressure drop over time when the anticoagulation level (heparin) is insufficient. Such a situation increases haemolysis rate and the likelihood of TMPa alarms. Carefully monitor the haemolysis and anticoagulation strategies in such circumstances.

18. Filter performance specifications require a minimum blood flow rate, specific to each filter, to avoid risk of hemoconcentration (see “Filter Operating Specifications” section).
19. During use, closely monitor the patient’s clotting parameters, especially when increasing the amount of anticoagulant delivered or after changing the anticoagulant syringe.
20. To assure proper anticoagulant flow control, only use the syringes listed in the operator’s manual. The use of non-recommended syringes can be a hazard for the patient. Particularly if there is no Luer-lock on the syringe, the seal between the syringe and the heparin line can no longer be guaranteed.
21. When not using the pre blood pump infusion circuit, it is recommended to clamp this circuit close to its connection to the access line; this will prevent the sedimentation of blood into the pre blood infusion line.
22. Always inspect the blood flowpath for signs of clotting before returning the blood in the set to the patient. If clotting is suspected, do not return the blood to the patient.
23. In case of re-circulation mode, the set must be replaced if the maximum re-circulation time is exceeded; refer to the PRISMAFLEX operator’s manual for more information.
In case of poor blood return, the set must be replaced.
In all cases, it is essential to reprime the set with fresh saline immediately before patient connection.
24. The PRISMAFLEX Set offers a specific design of the deaeration chamber which aims to trap air before blood is returned to the patient.
25. The PRISMAFLEX Set is not designed for a heater to be connected to the replacement solution line. A heater generates air bubbles which collect in the return deaeration chamber. Therefore, it is recommended not to use a heater on the replacement solution line. However, during treatment it may be necessary to warm the patient because of hypothermia.
26. The blood return line (blue-striped) is equipped with a Luer-lock connection near the deaeration chamber.
This connector is intended to join the extension line of a blood warmer. Refer to the specific Instructions for Use and strictly follow detailed instructions for set up of this line. Do not use this connection for any other purpose.
27. Do not attach/connect the extension line of a blood heater to the return line downstream of the air detector. The PRISMAFLEX system cannot detect air introduced in the line downstream of the air detector.
28. To minimize the risk of hemolysis, the PRISMAFLEX system monitors the TMPa and issues alarms if maximum pressure limits are exceeded. When performing TPE, additional monitoring for hemolysis is also recommended (see special procedures section on hemolysis).

 Warnings

1. **The use of operating procedures other than those published by the manufacturer or the use of accessory devices not recommended by the manufacturer can result in patient injury or death.**
2. Use only PRISMAFLEX Sets with the PRISMAFLEX system. The use of non-PRISMAFLEX Sets can result in patient injury or death.
3. Should acute allergic reactions (first-use syndrome) occur in patients receiving treatment via the PRISMAFLEX TPE1000/TPE2000 Set, **immediately stop the treatment** and administer appropriate intervention. Pay special attention to patients having already shown similar allergic reactions (see “Hypersensitivity Reactions” section).
4. Do not allow air to enter the blood compartment of the filter after priming is started. If a large amount of air enters, the set must be replaced.
5. Since drugs will pass through the membrane of the filter, the dosage of associated drug treatments must be adjusted for patients on TPE therapy.
6. To avoid complications due to increase of blood viscosity, the post-filler hematocrit (hematocrit in patient blood return) should never exceed 65%.
7. In some cases hemolysis may occur (pink or red tinge in the effluent bag) but may not be detected by a PRISMAFLEX alarm. When performing TPE, monitor the effluent bag closely (see “Special Procedures” section on “Hemolysis”).
8. Only use this set with the PRISMAFLEX system in TPE mode. Use of this set with any CRRT therapy can result in patient injury or death.

SPECIFICATIONS

See Tables at end of document.

SET MATERIALS

Plasmafilter hollow fiber	: Polypropylene
Housing and headers	: Polycarbonate
Potting compound	: Polyurethane
Tubing material	: Plasticized polyvinyl chloride (PVC)
Cartridge	: PETG

Note: the following information is available from the manufacturer upon request :

- information about test methods used to obtain performance characteristics,
- the number and range of particles in the effluent from the dialyzer prepared as recommended for clinical use,
- the types and amounts of residue from the sterilization process.

The PRISMAFLEX Set is not made with rubber natural latex.

Pressure monitoring in PRISMAFLEX TPE : definitions and purpose

*** access transmembrane pressure : TMPa = PA-PF**

TMPa is the local transmembrane pressure at the filter inlet and matches with the maximum local transmembrane pressure along the filter. It is used to predict risk of hemolysis, as hemolysis is associated with a threshold value in local filtration flux (<=> pressure gradient).

**** transmembrane pressure : TMP = (PA+PO)/2-PF**

This parameter represents the average pressure gradient across the membrane when linear blood and filtrate pressure drops are assumed. It allows to predict the filtrate flow rate.

(PA = filter inlet pressure (blood inlet side), PO = filter outlet pressure (blood return side), PF = filtrate pressure).

INSTRUCTIONS FOR USE

Note: use the set by following the detailed on-line instructions provided by the PRISMAFLEX system. Additional information is available in the PRISMAFLEX system operator’s manual.

Perform the following procedures when the appropriate instructions appear on the display of the PRISMAFLEX control unit.

		CAUTION	WARNING
LOAD SET	See online instructions		
PREPARE AND CONNECT SOLUTIONS	See online instructions	See caution n° 10	
PRIME SET	See online instructions Note: Priming solution: Saline or alkaline solution (pH>7.3) with added 5000IU heparin/liter.	See cautions n° 12 to 15 Cautions n° 26 and 28 Note : the PRISMAFLEX Set must be carefully deaerated	See warning n° 4

Anticoagulation Considerations

Note: see Cautions no. 18 through 21 and Warning no. 6.

Initiate anticoagulation of the blood flowpath, as prescribed by the physician. During use, monitor the patient’s clotting parameters; adjust the anticoagulation settings on the PRISMAFLEX system, according to the physician’s prescription. If prescribed, do not forget to infuse a loading dose of anticoagulant immediately after patient connection.

** Influence of Different Heparin Concentrations on the results of In Vitro Investigations in Plasma separation technology using Capillary Membrane Filters. J.K.Unger, C.Halterm... Artificial Organs 2003;27(7):649-647*

BAXTER CONFIDENTIAL - INTERNAL USE ONLY			
Part Number: 1000014865		Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.		Page: 8 of 24	
Colour Reference:		BLACK	

Anticoagulation plays an important part in extending filter life by retarding plugging and clotting.

		CAUTION	WARNING
CHANGE SET PROCEDURE	See online instructions	See caution n° 23	
RE-CIRCULATION PROCEDURE	See online instructions	See cautions n° 23 and 24	
END TREATMENT PROCEDURE	See online instructions	See caution n° 23	

MANUAL TERMINATION

Manual termination may be necessary due to power loss or an alarm of the PRISMAFLEX control unit. The alarm screen tells the operator if a manual termination is required.

Note: the following instructions are also found in “Troubleshooting” in the PRISMAFLEX system operator’s manual.

A. With Blood Return

Note: see Caution no. 23.

Note: a sterile spike connector may be required.

- Turn off the power. Clamp the access line (red-striped) and disconnect from the patient. Attach the access line to a 1-liter bag of sterile saline (use spike connector if needed). Unclamp the access line.
- Press the return clamp button (left side of the return line clamp assembly) and hold in the “in” position. With the other hand, remove the return line (blue-striped) from the return line clamp.
- Visually check the fluid level in the deaeration chamber. If the level is insufficient:
 - disconnect the deaeration chamber service line from the return pressure port on the PRISMAFLEX machine (the level will automatically rise in the deaeration chamber),
 - reconnect the line once the correct fluid level is reached.
- Remove the pump crank from its holder on the rear panel. Insert crank into the rotor of the blood pump and turn clockwise until sufficient blood is returned to the patient.

 **Warning:** the alarm system is disabled. Visually check for air in the blood return line until the patient is disconnected.

- Clamp the return line (blue-striped) and disconnect from the patient. Clamp lines to all bags.
- Press the two clips of the cartridge carrier to release the cartridge. Starting with the peristaltic pump, insert the pump crank into the rotor and turn each pump counterclockwise.
- When the pump segments are free, grasp the cartridge and pull out to disengage the lines from the pinch valves. Take the set off the control unit and discard as usual.

Note: remaining solutions may be used with a new set, if desired.

 **Warning:** ensure patient is disconnected from set before removing set from control unit.

B. Without Blood Return

Note: the patient will lose the blood contained in the blood flowpath during a manual termination without blood return.

- Turn off the power. Clamp the access line (red-striped) and return line (blue-striped) and disconnect from the patient.
- Clamp lines to all bags.
- Press the two clips of the cartridge carrier to release the cartridge. Starting with the blood pump, insert the pump crank into the rotor and turn each pump counterclockwise.
- When the pump segments are free, grasp the cartridge and pull out to disengage the lines from the pinch valves. Take the set off the control unit and discard as usual.

SPECIAL PROCEDURES IN CASE OF COMPLICATION

Filter Membrane Blood Leaks

Blood leaks through the filter membrane are automatically detected by the PRISMAFLEX control unit alarm system. A warning alarm is generated and blood loss is limited by immediate stoppage of all pumps.

To return blood to the patient, press STOP from the alarm screen, then press CHANGE SET from the Stop screen and follow the on-line instructions.

External Blood Leaks

Note: see Cautions no. 16, 17 and 21.

External blood leakage may not be immediately identified by monitoring equipment and could result in significant blood loss. Check the filter and all connections of the disposable tubings during treatment to minimize the risk of leakage. If an external blood leakage is observed, immediately stop the blood pump. Initiate corrective action by securing connections or replacing the PRISMAFLEX Set.

If necessary, administer adequate replacement solution to the patient to compensate for blood loss.

Hypersensitivity Reactions

Note: see Warning no. 3.

Should acute allergic reactions (first use syndrome) occur within the first few minutes of the treatment, it is important to react immediately by discontinuing the session and administering appropriate treatment.

Therefore special attention must be paid to patients having already shown similar reactions.

 **Warning:** blood in the extracorporeal circuit must not be returned to the patient.

Hemolysis

Note: see caution n°29 and Warnings n°6 and 7.

During TPE therapy, the risk of hemolysis exists due to the relatively large pore size of the membrane and the pressure difference needed between the blood and the effluent sides of the plasmafilter. The combination of a low blood flow rate and a high replacement flow rate (and or a high patient plasma loss flow rate) increases the risk of hemolysis by causing an excessive TMPa (access transmembrane pressure). To minimize the risk of hemolysis, the PRISMAFLEX system monitors the TMPa and issues alarms if maximum pressure limits are exceeded.

In some cases hemolysis may occur and not be detected by a PRISMAFLEX system alarm, however, a pink or red tinge may appear in the effluent bag. Therefore, while performing TPE therapy, monitor the effluent bag closely. In case of suspected hemolysis, immediately reduce the replacement and plasma loss flow rates to zero and verify patient’s parameters before proceeding with treatment.

WARRANTY AND LIMITATION OF LIABILITY

- The manufacturer warrants that the PRISMAFLEX Set has been manufactured in accordance with its specifications and in compliance with good manufacturing practices, other applicable industry standards and regulatory requirements.

If provided with the lot/serial number of the defective product, the manufacturer will, by replacement or credit, remedy manufacturing defects in the PRISMAFLEX Set becoming apparent before the expiration date.
- The warranty under paragraph a) above is in lieu of, and to the exclusion of, any other warranty, whether written or oral, express or implied, statutory or otherwise, and there are no warranties of merchantability or other warranties, which extend beyond those described in paragraph a) above. The remedy set out above for manufacturing defects is the sole remedy available to any person due to defects in the PRISMAFLEX Set and the manufacturer shall not be liable for any consequential or incidental loss, damage, injury or expense arising directly or indirectly from the use of the PRISMAFLEX Set, whether as a result of any defect therein or otherwise.
- The manufacturer shall not be liable for any misuse, improper handling, non-compliance with warnings and instructions, damage arising from events after the manufacturer’s release of the PRISMAFLEX Set, failure or omission to inspect the PRISMAFLEX Set before use in order to ensure that the PRISMAFLEX Set is in proper condition, or any warranty given by independent distributors or dealers.
- The manufacturer is GAMBRO Industries, 7 avenue Lionel Terray, BP 126, 69883 MEYZIEU CEDEX, FRANCE.

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 9 of 24	
Colour Reference:	BLACK	

INDIA

eng

- **Name and address of importer:** Baxter India Pvt. Ltd., Gala No. 1 to 6, Building No. C15, Ground Floor, Shree Arihant complex, Kalher, Bhiwandi, Dist. Thane, Tal: Bhiwandi (Thane-Zone 5) 421302
- Name, address, telephone no., email address (in case of consumer complaint):
Baxter India Pvt. Ltd.
5th Floor, Block A, Building 9,
DLF Phase III, Cyber City, Gurgaon- 122002.
Consumer Care No.: 0124-4603200
Consumer Care email ID: Customerservice_SHS_India@baxter.com
- Import License No.: MD-1591-2310

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 10 of 24	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 13 of 24	
Colour Reference:	BLACK	

GENERAL DATA /一般資料

	PRISMAFLEX TPE1000 SET	PRISMAFLEX TPE2000 SET
Weight / 重量	680 g	760 g
Blood volume in set / 套組中的血液容積 (± 10%)	71 mL	125 mL

FILTER OPERATING SPECIFICATIONS / 過濾器操作規格

	PRISMAFLEX TPE1000 SET	PRISMAFLEX TPE2000 SET
Minimum blood flow rate / 最低血液流速	50 mL/min	100 mL/min
Maximum blood flow rate / 最高血壓	180 mL/min	400 mL/min

	PRISMAFLEX TPE1000 SET	PRISMAFLEX TPE2000 SET	
TMPa / PTMa* maximum (mmHg) TMPa / PTMa* 最大值 (mmHg) QB/QS = 50 mL/min QB/QS = 100 mL/min QB/QS = 180 mL/min QB/QS = 200 mL/min QB/QS = 250 mL/min QB/QS = 400 mL/min	100 mmHg - 13,3 kPa 140 mmHg - 18,7 kPa 190 mmHg - 25,3 kPa –	– 120 mmHg - 16 kPa – 171mmHg - 22,8 kPa 193 mmHg - 25,7 kPa 246 mmHg - 32,8 kPa	
	180 mL/min	250 mL/min	400 mL/min
In vitro hemolysis rate (bovine blood, Htc 32%, Cp 60 g/L, 37°C)	< 1 mL blood/hour	< 1 mL blood/hour	< 3 mL blood/hour
體外溶血速率 (牛,血 , Htc 32%, Cp 60 g/L, 37°C)	< 1 mL 血液/h	< 1 mL 血液/h	< 3 mL 血液/h

FILTER DATA	過濾器資料	PRISMAFLEX TPE1000 SET	PRISMAFLEX TPE2000 SET
NOMINAL PHYSICAL CHARACTERISTICS Effective surface area Fiber internal diameter (wet) Fiber wall thickness	標稱物理特性 有效表面面積	0.15 m²	0.35 m²
	纖維內徑 (潮濕)	330 µm	
	纖維壁厚	150 µm	
IN VITRO PERFORMANCES ≠ Blood priming volume	體外效能 ≠ 血液預充量	23 mL ± 10 %	41 mL ± 10 %
Blood pressure drop (bovine blood, Htc****32%, Cp***** 60 g/L 37°C) QS/QB** = 50 mL/min QS/QB** = 100 mL/min QS/QB** = 180 mL/min QS/QB** = 200 mL/min QS/QB** = 250 mL/min QS/QB** = 400 mL/min	血壓下降 (牛血, Htc****32%, Cp***** 60 g/L, 37°C) QS/QB** = 50 mL/min QS/QB** = 100 mL/min QS/QB** = 180 mL/min QS/QB** = 200 mL/min QS/QB** = 250 mL/min QS/QB** = 400 mL/min	± 30 % 14 mmHg - 1,9 kPa 29 mmHg - 3,9 kPa 52 mmHg - 6,9 kPa – – –	± 30 % – 26 mmHg - 3,5 kPa – 51 mmHg - 6,8 kPa 64 mmHg - 8,5kPa 103 mmHg - 13,5kPa
Sieving coefficient (in vivo data - 19 treatments) - Albumin - IgG - Apolipoprotein B - IgM	篩濾係數 (體內數據 - 19 次治療) - 白蛋白 - IgG - 脂蛋白元 B - IgM	0.97 1 0.95 0.92	

≠ Typical mean values obtained from laboratory testing of post-sterilization sample lots. Results may vary depending on patient and clinical conditions. / 由實驗室測試滅菌後樣品所得之平均值。視病人及臨床狀況而定，結果可能不同。

* Transmembrane pressure / 跨膜壓 / ** Arterial blood flow rate / 動脈血液流速 / **** Hematocrit / 血容積比 / ***** Protein concentration / 蛋白浓度

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 14 of 24	
Colour Reference:	BLACK	

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BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 17 of 24	
Colour Reference:	BLACK	

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	PRISMAFLEX TPE1000 SET	PRISMAFLEX TPE2000 SET
Weight / 무게	680 g	760 g
Blood volume in set / 세트의 혈액량	71 mL	125 mL

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

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 18 of 24	
Colour Reference:	BLACK	

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
BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 19 of 24	
Colour Reference:	BLACK	

Prismatherm

SP420 line connexion



BAXTER CONFIDENTIAL - INTERNAL USE ONLY	
Part Number: 1000014865	Date: 06-JAN-20 Proofread No.: P03
Designer: E.M.	Page: 20 of 24
Colour Reference:	BLACK




INSTALL SET DIRECTLY ONTO

prismaflex

MONITOR

BAXTER CONFIDENTIAL - INTERNAL USE ONLY		
Part Number: 1000014865	Date: 06-JAN-20	Proofread No.: P03
Designer: E.M.	Page: 21 of 24	
Colour Reference:	BLACK	

Patented ready to use PrismaFlex set packaging.