

HEMODIALYSIS CONCENTRATES

Baxter is your one-stop
source for concentrates



Baxter

HEMODIALYSIS CONCENTRATES

Baxter Expands Portfolio: Baxter and Rockwell Medical



Baxter is now the exclusive distributor of Rockwell's hemodialysis concentrate products in the US. Rockwell is a leading manufacturer of products targeting chronic kidney disease with innovative products and services. This collaboration expands Baxter's product portfolio with the addition of Rockwell's hemodialysis concentrate products.

Citric acid concentrate powders and liquids

CitraPure acid concentrate is formulated to be used in conjunction with RenalPure Powder Bicarbonate (CB-Series) or SteriLyte Liquid Bicarbonate (CL-Series) only in a 45X three-stream hemodialysis machine. The listed electrolyte values represent the final dialysate electrolyte yield.

Formula	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Citric Acid (mEq/L)	Dextrose (mg%)
45X CitraPure Liquid Acid Concentrate - Drum (55-gallon)						
L410855	100.00	2.00	2.00	1.00	2.40	100.00
L411055	100.00	3.00	2.00	1.00	2.40	100.00
L411655	100.00	2.00	2.50	1.00	2.40	100.00
L411755	100.00	3.00	2.50	1.00	2.40	100.00
L412455	100.00	2.00	3.00	1.00	2.40	100.00
45X CitraPure Liquid Acid Concentrate - Case (4 - 1 gallon (3.78 L) containers/case)						
L41064	100.00	1.00	2.00	1.00	2.40	100.00
L41084	100.00	2.00	2.00	1.00	2.40	100.00
L41104	100.00	3.00	2.00	1.00	2.40	100.00
L41144	100.00	1.00	2.50	1.00	2.40	100.00
L41164	100.00	2.00	2.50	1.00	2.40	100.00
L41174	100.00	3.00	2.50	1.00	2.40	100.00
L41244	100.00	2.00	3.00	1.00	2.40	100.00
L41264	100.00	3.00	3.00	1.00	2.40	100.00
45X CitraPure Dry Acid Concentrate - Case (each case produces 25 gallons of CitraPure Liquid)						
D4108	100.00	2.00	2.00	1.00	2.40	100.00
D4110	100.00	3.00	2.00	1.00	2.40	100.00
D4116	100.00	2.00	2.50	1.00	2.40	100.00
D4117	100.00	3.00	2.50	1.00	2.40	100.00
D4124	100.00	2.00	3.00	1.00	2.40	100.00

Acetic acid concentrate powders and liquids

RenalPure and Dri-Sate acid concentrates are formulated to be used in conjunction with RenalPure Powder Bicarbonate (CB-Series) or SteriLyte Liquid Bicarbonate (CL-Series) only in a 45X three-stream hemodialysis machine. The listed electrolyte values represent the final dialysate electrolyte yield.

Formula	Na+ (mEq/L)	K+ (mEq/L)	Ca++ (mEq/L)	Mg++ (mEq/L)	Cl- (mEq/L)	Acetic Acid (mEq/L)	Dextrose (mg%)
45X RenalPure Liquid Acid Concentrate - Drum (55-gallon)							
C10855	100.00	2.00	2.00	1.00	105	4.00	100.00
C11055	100.00	3.00	2.00	1.00	106	4.00	100.00
C11455	100.00	1.00	2.50	1.00	104	4.00	100.00
C11655	100.00	2.00	2.50	1.00	105	4.00	100.00
C11755	100.00	3.00	2.50	1.00	106	4.00	100.00
C12455	100.00	2.00	3.00	1.00	106	4.00	100.00
C-13055	100.00	2.00	2.25	1.00	105	4.00	100.00
C-13155	100.00	3.00	2.25	1.00	106	4.00	100.00
45X RenalPure Liquid Acid Concentrate - Case (4 - 1 gallon (3.78 L) containers/case)							
C1024	100.00	0.00	0.00	1.00	101	4.00	100.00
C1064	100.00	1.00	2.00	1.00	104	4.00	100.00
C1084	100.00	2.00	2.00	1.00	105	4.00	100.00
C1104	100.00	3.00	2.00	1.00	106	4.00	100.00
C1144	100.00	1.00	2.50	1.00	104	4.00	100.00
C1164	100.00	2.00	2.50	1.00	105	4.00	100.00
C1174	100.00	3.00	2.50	1.00	106	4.00	100.00
C1244	100.00	2.00	3.00	1.00	106	4.00	100.00
C1264	100.00	3.00	3.00	1.00	106	4.00	100.00
45X Dri-Sate Dry Acid Concentrate - Case (each case produces 25 gallons of liquid acid)							
DC108	100.00	2.00	2.00	0.75	105	4.00	100.00
DC116	100.00	2.00	2.50	0.75	105	4.00	100.00
DC117	100.00	3.00	2.50	0.75	106	4.00	100.00

Bicarbonate concentrate powders and liquids

RenalPure Powder Bicarbonate and SteriLyte Liquid Bicarbonate are for use with CitraPure Dry or Liquid Acid Concentrate, Dri-Sate Dry Acid or RenalPure Liquid Acid Concentrate.

Formula	Na+ (mEq/L)	Bicarbonate (mEq/L)	Qty/case	Volume produced per bag (gallons)
RenalPure Powder Bicarbonate Bags				
CB0220C	37.00	37.00	20	2.1
CB154C	37.00	37.00	4	15.0
CB252C	37.00	37.00	2	25.0
SteriLyte Liquid Bicarbonate - Case (4 - 1 gallon (3.78 L) containers/case)				
CL024	37.00	37.00	4	NA

Custom formulations for 45X are available upon request. Other formulations for 35X and 36.83X are also available upon request. Contact your sales representative for more information.

CitraPure acid concentrate, Dri-Sate acid concentrate, RenalPure acid concentrate, RenalPure Powder

Bicarbonate and SteriLyte Liquid Bicarbonate are indicated for use in acute and chronic hemodialysis and to be used with the appropriate three stream hemodialysis machine/monitor and bicarbonate concentrate.

Rx Only. For the safe and proper use of the concentrates mentioned herein, refer to the Instructions for Use.

Rockwell provides the following White Glove Services so that customers receive the special care and products they deserve.



Dry Acid Concentrate Mix System-DAMX45

- Mixes 100 gallons of product in approximately 3 minutes
- Tank = 34" in diameter
- Distance from base to RO (reverse osmosis) Fill Port (where product is added) = 48" high
- Uses forced RO water pump to mix the product



For more information visit www.esrddialysis.com,
or contact your Baxter Sales Representative



Customer Service*

- Tracks all orders and deliveries
- Ensures stock levels are adequate in cases of non-regularly scheduled orders
- Verifies all codes to prevent ordering errors
- Selects most cost effective delivery option for urgent orders (common carrier, drop ship, etc.)
- Calculates minimum order requirements to avoid extra fees

*Ranked #3 on the DaVita annual supplier listing



Transportation

- Pump-off services for drum users
- Drivers remove and recycle empty drums
- Deliver product to customer's desired storage location
- Target date and time for delivery is provided
- Transportation fleet is equipped with satellite tracking equipment



Products

- RENALPURE liquid acid jugs are 3.78 L vs. FMC NATURALYTE liquid acid jugs which are 3.43 L
- CITRAPURE is the only 100% acetate-free product in the U.S. market
- RENALPURE bicarbonate (CB0220C) contains 20 bags per case vs. FMC NATURALYTE bicarbonate (O2-4400-1) containing 12 bags
- RENALPURE bicarbonate (CB252C) contains two 25 gallons bags vs. FMC NATURALYTE bicarbonate containing 1 bag